



Streamlining Ostomy Care

Designed in Canada.

Manufactured in China for OstoVersa Care Ltd. British Columbia, Canada.

LiberOsto™ Cleaning and Disinfection

Cleaning and Disinfection

1. Rinse the Funnel and Pipe immediately after use with clean water.
2. If necessary, wash the Funnel and Pipe with a mild liquid soap and soft toilet paper, then discard the used toilet paper into the toilet.
3. Disinfect the Funnel and Pipe before overnight storage after use; after each use during an infectious illness; after transport if contamination is suspected; and before long-term storage.
4. Use either 0.1% sodium hypochlorite (bleach), prepared according to the manufacturer's instructions and applied for 10 minutes, or a quaternary ammonium disinfectant used as directed by its manufacturer.
5. Apply disinfectant to Funnel and Pipe surfaces exposed to ostomy output or other contamination. Rinse when required by the disinfectant instructions.
6. Allow the Funnel and Pipe to drain and air-dry completely.

7. If contaminated, wipe or rinse the Shaft Assembly, Pivoting Bracket, and Base Assembly with clean water as needed. Do not immerse these components in a bathtub.
8. Do not replace the Cover until the Funnel is completely dry.
9. Do not use water hotter than your hands can comfortably tolerate (typically no more than 40 °C / 104 °F), abrasive cleaners, abrasive brushes or scouring pads, sharp objects, or strong solvents.

12. Disinfection warnings:

1. Do not mix bleach with ammonia, vinegar, acids, alcohol, or other cleaning products.
2. Use disinfectants in a well-ventilated area, wear any protective equipment recommended by the disinfectant manufacturer.
3. Do not use disinfectants other than those specified in the Detailed IFU.

13. Other Warnings

1. ⚠ Do not use LiberOsto™ if any component is missing, defective, or damaged.
2. ⚠ Stop using the device immediately if its height or angle changes unexpectedly, or if any connection becomes loose during use.
3. ⚠ Do not lean on the device or place your body weight on it.
4. ⚠ Ensure that the Funnel and Pipe are completely empty before adjusting or releasing the angle-locking mechanism.
5. ⚠ Keep fingers away from the junction between the moving Pivoting Bracket and the Shaft Assembly during angle adjustment.
6. ⚠ Keep the assembled device, detached components, packaging materials, and all cleaning and disinfecting products out of the reach of children.